

IS YOUR LAB READY TO REV ITS ENGINE?

Here Are a Few Reminders For a Successful Return to Cell Cultivation

You wouldn't start a car that has been stationary for months without a maintenance check. Don't start waking up those cells without completing these steps!




STEP 1 PREPARE TO INITIATE



- Make sure to turn off the power to your incubator and the humidity pan is empty.
- Take the time to complete a full interior wipe down with 70% EtOH.



STEP 2 DECONTAMINATE




- Run the appropriate decontamination cycle for your incubator:

Vaporized Hydrogen Peroxide Cycle, 180°C High Heat Sterilization Cycle, 24 hour UV Cycle.



STEP 3 INVESTIGATE AND CALIBRATE



- Run the incubator at 37°C for 24 hours without gas.
- Check the gas levels of your CO₂ or N₂ tanks. Replace as needed.
- That gas won't go far if your tubes are frayed or cracked — so inspect and replace. And don't forget your inline filters too!
- You are now ready to open the gas supply.
- Run a gas calibration — or better yet, contact a service professional to complete the procedure to ensure calibrated gas levels.



Now you are set to wake up your cells and ESBE is here to help!

Ready for a new incubator? Ready to culture in 3D?

Call your local ESBE representative to learn more about our portfolio of Cell-IQ™ incubators and PrimeSurface® culturing well plates.

CELL CULTURING
JUST GOT BETTER

CELL-IQ™ INCUBATORS



Cell-IQ™ Series

PAN-MCO170AICUVHLPA

5.8 cu.ft.

CO2 Laboratory Incubator



Cell-IQ™ Series

PAN-MCO230AICUVHLPA

8.1 cu.ft.

CO2 Incubator for Cell Culture



Cell-IQ™ Series

PAN-MCO170MUVHPA

5.7 cu.ft.

Multigas CO2/O2 Incubator

BACK TO THE LAB PROMOTION

Get free Atomizer and Roller Base
when you purchase incubator(s).

Terms & Conditions Apply.
Promo valid until March 31, 2021.

Contact your local ESBE Representative
for more information.

FREE Items



PAN-MCOHPPW(PROMO)
H2 O2 Generator/ Atomizer, \$1575 value
FREE with the purchase of one of the above listed incubators.
Limited to one free atomizer per lab.



PAN-MCO170RBPW(PROMO)
Roller base for Cell IQ Touch (KM-CC17RU2A), \$1030 value
FREE with the purchase of two stacked incubators; qualifying models:
PAN-MCO170AICUVHLPA and PAN-MCO170MUVHPA.



PAN-MCO230RBPW(PROMO)
Roller Base for MCO-230AIC, \$1260 value
FREE with the purchase of two stacked incubators;
qualifying model: PAN-MCO230AICUVHLPA.