

 **OWEN MUMFORD**
Making a World of Difference

www.owenmumford.com

UK / EU

Owen Mumford Ltd.
Brook Hill, Woodstock
Oxford OX20 1TU, United Kingdom
T : 01993 812021
F : 01993 813466
E : info@owenmumford.net

NORTH AMERICA

Owen Mumford Inc.
1755 West Oak Commons Court
Marietta, Georgia 30062, USA
T : 800-421-6936
F : 770-977-2866
E : info@owenmumfordinc.com

FRANCE


Owen Mumford
BP 444, 27204 Vernon Cedex, France
T : 02 32 51 88 70
F : 02 32 21 96 95
E : contact@owenmumford-fr.com

GERMANY

Owen Mumford GmbH.
Alte Häge 1
63762 Großostheim-Ringheim, Germany
T : 06026/97750
F : 06026/977529
E : info@owen-mumford.de

EXPORT

Owen Mumford Ltd,
Brook Hill, Woodstock
Oxford OX20 1TU, United Kingdom
T : +44 1993 812021
F : +44 1993 813473
E : export@owenmumford.co.uk



Autoject® 2 is a registered trademark of Owen Mumford Limited.
All other trademarks are the property of their respective owners.



Autoject® 2



The simple solution



AJ/1300/MCC/0607/7

Autoject[®] 2

Self-injection made easy

Autoject 2 automates the process of self-injection. At the touch of a button Autoject 2 injects the needle to a pre-set depth and delivers the contents of the syringe.

There are two standard Autoject 2 models:

1. For patients using fixed needle syringes.
2. For patients using removable needle syringes.

Each Autoject 2 comes in a handy zipped soft pouch that includes:

- 1 Autoject 2
- 2 depth adjusters (FIXED NEEDLE)
- 1 depth adjuster (REMOVABLE NEEDLE)
- a yellow primer/depth guide
- instructions for use




Fixed Needle Autoject 2 set shown above

Raised trigger button for firing
Easy for both children and seniors to operate.

Window showing the injection is completed
When the trigger button is pressed, patients know the dose has been fully delivered by checking for the yellow marker that can be seen through the window on the injector body, giving the assurance that the injection has been completed.

Safety interlock mechanism
When the Autoject is loaded with a syringe, the safety device prevents accidental firing until the device is pressed against the skin.

A choice of depth adjusters
Allows the use of a wide range of syringe brands and the ability to pre-set different penetration depths.



COMFORT

Nobody likes injections, so for people who have to inject, comfort is vital. The Autoject 2 range offers easy to handle devices for quick injections to reduce discomfort.

CONFIDENCE

Many people choose Autoject 2 because they are concerned by injections, and particularly injecting themselves. Autoject 2 offers a perfect solution by concealing the needle from the patient's view prior to injection.

Autoject[®] 2 Fixed needle

AJ1300



Compatible with most major brands of fixed needle syringe.

Autoject[®] 2 Removable needle


AJ1311



Specifically designed for compatibility with the 1ml BD Tuberculin syringe.

Autoject[®] EI Fixed needle

AJ1310



Compatible with most brands of fixed needle syringes. Autoject EI's External Indicator provides confirmation that the full dose has been delivered.