


# Cellometer<sup>®</sup> Auto 1000

## Cell Counter

### Accurate Cellometer Cell Counting ... Now Even Faster and Easier

- Improve counting of live/dead cells, individual cells within clusters, and irregular-shaped cells: Proprietary Cellometer software
- Spend less time counting and more time on research: Concentration and viability in <10 seconds!
- Easily count cells in the hood: Stand-alone system with a large, user-friendly touch screen
- View, print, and save counted cell images, data reports, and cell size histograms: Connect to a network or use a USB drive




## Simple, Automated Cell Counting


Improve Speed, Accuracy, and Consistency



### Load 20 $\mu$ l of Sample




### Insert Counting Chamber



### Select Assay & Click Count



### View Image & Results



For more information on the new Cellometer Auto 1000, contact Nexcelom at 978-327-5340 or [info@nexcelom.com](mailto:info@nexcelom.com)!