

# ROTINA 380 BENCHTOP CENTRIFUGE

FOR **CLINICAL & RESEARCH** APPLICATIONS

## FEATURES:

- Robust & steady, temperature control system with pre-cool function from -20 to +40 °C (Optional)
- Versa-Rotor™ for multiple applications (no need to exchange rotors)
- Run-Silent™ Operation ( $\leq 65$  dB, as soft as normal conversation)
- Easy-Lift™ / Safety-Lock™ Lid
- Five (5) Year Warranty
- FDA registered centrifuge manufacturer

**Hettich**

CENTRIFUGES

## CENTRIFUGE PACKAGES

**NEW**

- BLOOD
- CELL CULTURE
- PLATES
- BOTTLES

(See reverse for more info)



Find out more by visiting:  
[www.hettweb.com/rotina-380](http://www.hettweb.com/rotina-380)

MAX. CAPACITY:  
**4 x 290 mL**


**DIMENSIONS (HxWxD)** ROTINA 380: 16.4" x 18" x 23.6" | ROTINA 380R: 16.4" x 18" x 29.5"

**TUBE TYPES** MICROLITER, ROUND BOTTOM, BLOOD, URINE, KOVA, CONICAL, SKIRTED, SCREW TOP, BOTTLES, LARGE CONICAL, PLATES, ZENTRIFORM


# ROTINA 380 CENTRIFUGE PACKAGES FOR CLINICAL & RESEARCH APPLICATIONS

(BELOW IS A LIST OF THE **MOST POPULAR** CENTRIFUGE PACKAGES AVAILABLE FOR THIS UNIT. **ADDITIONAL CONFIGURATIONS** ARE AVAILABLE DEPENDING ON YOUR APPLICATION. PLEASE INQUIRE FOR MORE DETAILS)


## ROTINA 380 BLOOD TUBE PACKAGE - (4 or 2)

	Tube Qty	Volume (mL)	Size (mm)	RPM	G-FORCE	<b>PACKAGE CONTENTS</b>  (1) Centrifuge (1) 90° Swing-out Rotor (4-place) (4) Buckets (4) Adapters for 13 x 75/100 mm (4) Adapters for 16 x 75/100/125 mm <b>or</b> (2) Adapters for 13 x 75/100 mm (2) Adapters for 16 x 75/100/125 mm
	64 (32)	1.6 - 7	13 x 75/ 100	5,000	3,913	
	52 (26)	4 - 10	16 x 75/ 100			
	12 (6)	8	16 x 125			
<b>Pkg No.:</b>				<b>380BLOOD (4 or 2)</b>		
<b>Pkg No. w/Temperature Control:</b>				<b>380RBLOOD (4 or 2)</b>		
<b>Add Bio-containment Lids</b>				<b>(-BC)</b>		


## ROTINA 380 CELL CULTURE PACKAGE - (4 or 2)

	Tube Qty	Volume (mL)	Size (mm)	RPM	G-FORCE	<b>PACKAGE CONTENTS</b>  (1) Centrifuge (1) 90° Swing-out Rotor (4-place) (4) Buckets (4) Adapters for 15 mL (4) Adapters for 50 mL <b>or</b> (2) Adapters for 15 mL (2) Adapters for 50 mL
	36 (18)	15   Conical	17 x 120	5,000	4,863	
	16 (8)	50   Conical	29 x 115			
	<b>Pkg No.:</b>				<b>380CELLCULTURE (4 or 2)</b>	
<b>Pkg No. w/Temp. Control:</b>				<b>380RCCELLCULTURE (4 or 2)</b>		
<b>Add Bio-containment Lids</b>				<b>(-BC)</b>		

## ROTINA 380 PLATE PACKAGE

	Tube Qty	Volume (mL)	Size (mm)	RPM	G-FORCE	<b>PACKAGE CONTENTS</b>  1) Centrifuge (1) 90° Swing-out Rotor (4-place) (2) Buckets
	10	Microtiter	--	4,000	2,397	
	<b>Pkg No.:</b>				<b>380PLATE</b>	
<b>Pkg No. w/Temperature Control:</b>				<b>380RPLATE</b>		

## ROTINA 380 BOTTLE PACKAGE

	Tube Qty	Volume (mL)	Size (mm)	RPM	G-FORCE	<b>PACKAGE CONTENTS</b>  1) Centrifuge (1) 90° Swing-out Rotor (4-place) (4) Buckets (4) Adapters
	4	250	61 x 122	5,000	4,863	
	<b>Pkg No.:</b>				<b>380BOTTLE</b>	
<b>Pkg No. w/Temperature Control:</b>				<b>380RBOTTLE</b>		

↓  
Find out more by visiting:  
[www.hettweb.com/rotina-380](http://www.hettweb.com/rotina-380)

DISTRIBUTED BY:

**Hettich**  
CENTRIFUGES

**HETTICH NORTH AMERICA**  
Toll Free: 1-866-370-4388  
Email: [info@hettweb.com](mailto:info@hettweb.com)

→ Discover more solutions online at: [www.hettweb.com/centrifuges](http://www.hettweb.com/centrifuges)

Copyright © 2012 Hettich Instruments, LP. All rights reserved. For the most up to date information on Hettich please visit [hettweb.com](http://hettweb.com).  
Printed in the USA using soy based inks on FSC certified recycled paper (25% pre-consumer waste & 25% post consumer waste).

PS-ROTINA-380-E RevA 01-2013