



Benchtop Clinical Centrifugation Solutions

StatSpin®

Now Available Through ESBE



ESBE Scientific

The StatSpin® Lineup



Express 2 Centrifuge

High-quality serum or plasma samples in only two minutes

ESBE Item ID: HEM-STSSX2



Express 3 Centrifuge

Achieve fast turnaround time and greater productivity

ESBE Item ID: HEM-STSSX3



Express 4 Centrifuge

Achieve gel separation in only three minutes

ESBE Item ID: HEM-STSSH4



Cytofuge 2 Cyto centrifuge

A quiet and compact cyto centrifuge system

ESBE Item ID: HEM-STCF02



MP Centrifuge

Blood, urine, hematocrits and other routine spins in two minutes or less

ESBE Item ID: HEM-STSSMP

Quality Sample Separations to Meet Growing Demand

The StatSpin family of high-speed tabletop centrifuges provides quality sample separation to meet the demand for rapid test turnaround time, delivering operational efficiency and staff satisfaction for your laboratory at a lower cost.

Experience rapid turnaround time and workflow efficiencies



Deliver high g-forces and provide two- to five-minute STAT sample processing for gel or non-gel tubes



Enhance workflow to achieve faster results and improve turnaround times for Emergency Department and other critical samples

Enhance staff satisfaction



Process samples on quiet tabletop centrifuges with simple operation and minimal maintenance requirements



Free your staff from spending unproductive time waiting for sample separation on a conventional centrifuge

ESBE Service



We provide full service offering for the Hemocue StatSpin® Line of Clinical Benchtop Centrifugation Solutions



Full Suite of In-house and On-site Services. Licensed and ISO certified service provider.



**StatSpin Express 2
Centrifuge**

High Productivity. Lower Cost.

The StatSpin Express 2 primary tube centrifuge prepares high-quality plasma or serum samples in only two minutes. Meet rapid test turnaround requirements without a major investment in equipment or a disruption in your laboratory routine.

ORDERING INFORMATION:

StatSpin Express 2 Centrifuge for 100-240 V, 50/60 Hz
Includes 4-place rotor (RTX4), universal power supply, removable chamber liner, four 5 mL inserts and North American power cord.

ESBE Item ID: HEM-STSSX2



[View Brochure](#)



**StatSpin Express 3
Centrifuge**

Faster turnaround times and greater productivity.

Place it next to any of your STAT analyzers and serum or plasma is ready in just minutes. The StatSpin Express 3 primary tube centrifuge can spin eight tubes up to 10mL. The centrifuge provides high-quality plasma or serum in two to five minutes.

ORDERING INFORMATION:

StatSpin Express 3 Centrifuge for 100-240 V, 50/60 Hz.
Includes 8-place rotor, removable bowl liner, eight 5 mL inserts and power cord.

ESBE Item ID: HEM-STSSX3



[View Brochure](#)



**StatSpin Express 4
Centrifuge**

Achieve gel separation in only three minutes.

The StatSpin Express 4 high-speed horizontal centrifuge provides fast turnaround time and greater productivity for your Lean laboratory. Place it next to any STAT analyzer, and serum or plasma gel tubes are ready in just minutes. The horizontal spin provides a flat gel barrier.

ORDERING INFORMATION:

StatSpin Express 4 Centrifuge for 100-240 V, 50/60 Hz
Includes 8-place rotor, eight 7 mL inserts, eight 10 mL inserts, manual latch release tool, one hex key tool, and power cord.

ESBE Item ID: HEM-STSSH4



[View Brochure](#)



StatSpin® Now Available through ESBE Scientific

Promotional pricing
available for a limited time



StatSpin Cytofuge 2
Cyto centrifuge

A quiet and compact cyto centrifuge system.

The StatSpin Cytofuge 2 is a personal cyto centrifuge that costs thousands of dollars less but provides the same high-quality sample preparations as more expensive systems. Its simple, understandable controls and snap-seal filter concentrators make operation easier than ever before.

ORDERING INFORMATION:

StatSpin Express 2 Centrifuge for 100-240 V, 50/60 Hz
Includes 4-place rotor (RTX4), universal power supply, removable chamber liner, four 5 mL inserts and North American power cord.

ESBE Item ID: HEM-STCF02



StatSpin MP
Centrifuge

Flexible solutions. Ready in two minutes.

The StatSpin MP centrifuge is used for centrifuging microhematocrit tubes for packed cell volume determination, as well as rapid separation of whole blood and preparation of urine sediment for microscopic analysis. Quality serum or plasma separations takes as little as 30 seconds and no more than two minutes.

ORDERING INFORMATION:

StatSpin Express 2 Centrifuge for 100-240 V, 50/60 Hz
Includes 4-place rotor (RTX4), universal power supply, removable chamber liner, four 5 mL inserts and North American power cord.



ESBE Item ID: HEM-STSSMP
with HCT reader

HEM-STB71954
without HCT reader